




TRANSPARENT FITTINGS



Precision manufactured, transparent fittings allow for user observation of critical junctions. These fittings offer a window to notoriously problematic microfluidic connections, providing researchers a better understanding of their experiment. Fittings come in a variety of architectures, port styles, and inner diameters. Don't see what you need? Please contact us at info@myCorSolutions.com.

Part Number	Architectures	Port Style	Inner Diameter
U-1/4-28-1.0	Union	¼-28 Flat Bottom	1.0 mm
U-1/4-28-1.6	Union	¼-28 Flat Bottom	1.6 mm
U-10-32-0.8	Union	10-32	0.8 mm
Y-1/4-28-1.0	"Y"	¼-28 Flat Bottom	1.0 mm
Y-1/4-28-1.6	"Y"	¼-28 Flat Bottom	1.6 mm
Y-10-32-0.8	"Y"	10-32	0.8 mm
T-1/4-28-1.0	"Tee"	¼-28 Flat Bottom	1.0 mm
T-1/4-28-1.6	"Tee"	¼-28 Flat Bottom	1.6 mm
T-10-32-0.8	"Tee"	10-32	0.8 mm
X-1/4-28-1.0	Cross	¼-28 Flat Bottom	1.0 mm
X-1/4-28-1.6	Cross	¼-28 Flat Bottom	1.6 mm
X-10-32-0.8	Cross	10-32	0.8 mm